Calcium-Trenn-Set

Aufgabe

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Calcium Separation Set

Purpose

The calcium separation set makes it possible to perform a cadmium determination with the Dr. Lange cuvette test LCK 308. Calcium and magnesium ions are selectively separated from cadmium ions.

Range of application

The separation set can be used to treat samples containing calcium and magnesium concentrations of up to 5000 mg/l so that interference-free cadmium determinations can be carried out. If higher concentrations are present the whole separation process must be repeated.

Contents of package

1. Molecular sieve
2. Membrane filtration attachments
3. Disposable syringe
4. Dosing spoon

Accessories

Dr. Lange cuvette test LCK 308
Beaker (HBG 010)

Storage Information

The separation set can be stored at +15 to 25°C for an indefinite period.

Information

The separation of calcium and magnesium is only possible in conjunction with the complexing solution A of the LCK 308 cuvette test. Other applications for the separation set have not yet been tested.

For special attention

After each separation the syringe must be thoroughly rinsed with distilled water.

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1. Separation

To a beaker add

pretreated water sample from LCK 308	complete contents of the reaction tube
Molecular sieve	1 spoonful

Shake beaker for 1 min, then allow to stand for at least 10 min.

2. Filtration

Draw the supernatant liquid into the syringe. Fit the filtration attachment firmly on the Luer connection of the syringe. Hold the syringe with the filtration attachment pointing upwards and slowly press the liquid upwards until the membrane filter is wetted. Press the sample slowly through the filter attachment. Discard the first 15 drops (approx. 1 ml) and catch the rest in the reaction tube from LCK 308.

The clear solution obtained is virtually free from calcium and magnesium ions and can be analysed for cadmium. The further analysis is carried out in accordance with the working procedure for cuvette test LCK 308 (see item 3. Analysis).

Dear user!

These instructions contain the analyses for all Dr. Lange equipment.

Please cut out the analysis for your equipment with scissors and file it in the Analysis File (Cat. no. BDZ 014). In this way you will be able to collect analyses which apply to your equipment.

Caution

The analyses which you collect in this way must always be updated whenever Dr. Lange changes the edition date on the instructions.

If the edition date on the enclosed slip does not agree with the edition date on the filed analysis, this means the filed analysis is no longer current.

Please replace it by the new analysis on the enclosed slip.

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